

The Moffatiser™


This product is for routine applications of local anaesthetic on the upper airways mucosa including anterior and posterior nasal space, oral cavity and the naso-lacrimal duct (DCR). The Moffatiser™ allows the physician to apply local anaesthetic in an efficient, accurate and sterile way. It is safe and easy to use.

Product features:

- Ergonomic design for comfortable and easy holding
- Malleable distal tip to allow accurate placement of local anaesthetic
- Sterile and disposable to eliminate the transfer of harmful pathogens between patients
- Smooth action syringe for gentle application

The Moffatiser™ is particularly useful for the following applications:

- ENT outpatient procedures
- Respiratory procedures
- Awake fibre optic intubation
- Upper ENT procedures



Code	Product Description	Pack Quantity
6100	Moffatiser™ + 5mL syringe + 18g drawing needle	20
6101	Moffatiser™ + 1.8m oxygen tubing	20